BEGO GUIDE PLANNING CENTER





Registration in the BEGO Guide order portal is naturally free of charge.

Do you have any questions? The team BEGO Guided surgery will happily help you with your inquiries under: +49 421 2028-230



EASY. SAFE. IMPLANTATION.



BEGO Implant Systems GmbH & Co. KG

Wilhelm-Herbst-Str. 1 · 28359 Bremen · Germany Tel. : +49 421 2028-230 · Fax : +49 421 2028-265 guide@bego.com · guide.bego.com

Always want to have a look at all BEGO news? Sign up for our newsletter here: www.bego.com/newsletter













BEGO Guide

Quick and easy fully guided surgery



84618/01 · IHa· © 2023 by BEGO · 20230530

GUIDED SURGERY — QUICK AND EASY

YOUR BENEFITS AT A GLANCE

PROCEDURE

BEGO Guide surgical guides are the optimal supplement to the BEGO Semados® implant system.

Maximum flexibility – low costs

- Pilot drill hole or fully guided
- Tooth or mucosa supported
- Optional transgingival, minimally invasive surgery possible
- Including digital wax-up on request*
- Upon request, design and fabrication of a temporary restoration based on the implant planning.*

No need to have your own planning software

Order your proposed plans in the BEGO Guide order portal – we'll do the planning for you!



Planning for implant systems from all manufacturers

Planning services for BEGO Semados® implants and implant systems from other manufacturers.

Pilot guides

- For the pilot drill hole
- Compatible with every implant system
- Sleeves with 2 mm diameter

Fully guided surgery

- For fully guided insertion of BEGO Semados® S/SC/SCX/RS/RSX implants
- Fitted with BEGO Guide master sleeves
- Guided drilling of all drilling steps
- Defined depth stop

Priceless.

BEGO Security Plus for even greater security

5 year warranty for restorations with BEGO Semados® implants and their prosthetic components

 Can be selected in the ordering process in the BEGO Guide order portal



Three steps to a BEGO Guide surgical guide:



Register in the BEGO Guide order portal

Register at guide.bego.com and you will receive your personal access data.



Create a patient case and upload the CT/CBCT data

- Create the case and provide details of the individual planning requirements
- Upload the CT/CBCT data to the portal or send in a data CD
- Ship the model using the collection service
- Your individual proposed plan is prepared



Expert review and sign-off of the plan

- Proposed plan is ready for viewing and expert review
- Online session possible for required changes
- Plan is signed off in the portal
- Guides are produced in Bremen
- Surgical guides are shipped in just a few days



^{*} optional scan and design service (additional charge)