



### Milling blank for the generation of occlusal splints with thermomemory effect

Milling blanks made of BEGO PMMA Splint E are characterized by a thermoplastic flexibility with thermal memory effect. Due to the industrial production process the highest material homogeneity is achieved, which guarantees outstanding long-term stability. The use of CAD/CAM technology also ensures a safe process because mixing errors (e. g. by hand mixing) are eliminated.

The special material properties result in highly precise adaptation to the tooth situation and extraordinary, stress-free wearing comfort for the patient. Furthermore, the self-adjusting bite splint is extremely fracture-proof and has a high optical transparency.



Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.

HIGH AND PRECISE

# ADAPTATION TO THE TOOTH SITUATION

### Special features at a glance and technical details

### The advantages for you

- High and precise adaptation to the dental bite situation
- Exceptional, tension-free functional comfort for the patient by thermomemory effect
- Extremely fracture-proof

- Highest thermomemory effect on market (returns to original shape between uses)
- High optical transparency
- No adverse taste

### **Product details**

### **Chemical composition**

Poly(m)ethylacrylate and cross-linking copolymers of methacrylic acid	> 90 %
1.2-cyclohexane dicarboxylic acid diisononyl ester	< 10 %

#### Material data

Flexural strength (23 °C)	> 20 MPa
Flexural strength (37 °C)	< 20 MPa
Self-alignment (37 °C)	> 95 %
Density	approx. 1.1 to 1.2 g/cm <sup>3</sup>
Color	transparent

Availability	Diameter	Contents	REF
Milling blank PMMA Splint E 20 mm	98.5 mm	1 Piece	71200
Milling blank PMMA Splint E 16 mm	98.5 mm	1 Piece	71201

**EXCEPTIONAL, TENSION-FREE FUNCTIONAL** 

# COMFORT

FOR THE PATIENT BY THERMOMEMORY EFFECT

### **Product description**

- The high-tech plastic is included in the group of smart materials\* and is designed to generate functionally therapeutic occlusal splints
- BEGO PMMA Splint E becomes flexible and thus self-adjusting from 27 °C thanks to its thermoactive property
- The material's thermomemory effect translates into the highest self-aligning torque at room temperature available on the market
- Splints made of BEGO PMMA Splint E are extremely tear- and fracture-proof, in addition to being long-lasting and permanently transparent
- The material's composition satisfies high medical standards. The BEGO PMMA Splint E does not contain any phthalate-based softeners; it only features flexibilizers approved for medical technology applications. This significantly reduces the material's irritation and sensitization potential
- The use of high-quality raw materials ensures a high level of biocompatibility and very good tolerability
- The unique production process eliminates any adverse effects on taste sensation
- BEGO PMMA Splint E satisfies the requirements of a Class IIa medical device as defined by Directive 93/42/FFC

<sup>\*</sup>Smart materials (also known as responsive materials) is a term used to refer in a narrower sense to materials that were specifically developed to independently respond to changing ambient conditions (such as temperature) in a specific way. This change is reversible.

# SELF-ADJUSTING, EXTREMELY FRACTURE-PROOF, HIGH OPTICAL TRANSPARENCY

### Milled splints by BEGO Medical

We also offer you milled splints from our milling center in Bremen. Send us your file of the designed splint (in \*.STL, \*.BEGOSTL, \*.3ox or \*.3oxz format) and we'll mill the splint in our very own production facilities. BEGO Medical customer support will be more than happy to provide you with the necessary design parameters.

### **General information**

User Support (for technical enquiries)		Contact and service (ordering service, invoicing queries)	
Service hours Mon. – Thu. Fri.	8:00 am – 6:00 pm 8:00 am – 5:00 pm	Service hours Mon.—Thu. Fri.	8:00am-5:00pm 8:00am-4:00pm
Contact Telephone E-Mail	+49 421 2028-200 info@bego-medical.de	Contact Telephone Fax E-Mail	+49 421 2028-220 +49 421 2028-240 sales@bego.com



Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.





### BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG

Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany Tel. +49 421 2028-0 · Fax +49 421 2028-100 E-Mail info@bego.com · www.bego.com



Pro3dure medical GmbH Am Burgberg 13 58642 Iserlohn

Always want to have a look at all BEGO news? Sign up for our newsletter here: www.bego.com/newsletter









